



## Connect Access to a SurveyMonkey survey

### 1. Create your survey

---

1. Create your survey in SurveyMonkey.
2. Click the tab **Collect Responses**.
3. Choose **Get Web Link**.
4. Copy the Web Link

### 2. Add Web Link to Access

---

1. Specify your project and target groups in Access.
2. Click **Create Project**
3. For each of your target groups, click on **Add links**.
4. **Live link** (*Will be used by your study in field*)
  - a. Paste the **web link** link from SurveyMonkey into **link template**.  
At the end of the link, add **?c=[ID]**. The link should look like this:  
[https://www.surveymonkey.com/XXX/XXXXXXXX?c=\[ID\]](https://www.surveymonkey.com/XXX/XXXXXXXX?c=[ID])
5. **Test link** (*Will be used by Cint to test your study before approving it*)
  - a. Paste the live link from step 4 into the **link template** field.
6. Enter **Screenout and testing instructions** for Cint
  - a. Do you use any screenout questions? If so, detail them along with the answers that lead to a screenout
7. Confirm your target group setup

### 3. Add redirect links to SurveyMonkey

---

1. In **SurveyMonkey**, go to the **Collect Responses** tab
2. Choose a **collector**
3. Click **Show advanced options**.
4. Make sure **Custom Thank you** is OFF
5. Select **Survey End Page** (mandatory)
  - a. Choose **On, show a custom end page upon survey completion** and enter the **complete link** from **Access**
6. Select **Custom Disqualification** (if your survey contains screenout questions)
  - a. Choose **On, send disqualified respondents to a custom disqualification URL** and paste the **screenout link** from **Access**
7. The settings in SurveyMonkey will be automatically saved

### Advanced – Screenout logic in SurveyMonkey

---

To use screen-outs in SurveyMonkey, use Question Logic and choose to skip to the Disqualification Page. This will screen out all respondents who have selected the "wrong" answer.