


How to set up your Survey Lab survey with Cint


SurveyLab works with different online survey panels like CINT, mypanellab.com, luc.id, or Dynata (ResearchNow/SSI).

In this article, you will find short instruction that describes how to conduct a survey using CINT survey panel.

Adding Cint's redirects

1. Log into SurveyLab
2. Create a survey collector (**URL Link**) or edit the existing one



3. Go to survey collector  **SETTINGS** and add redirection links (copied from CINT).

New collector

 **SETTINGS**  **SHARE**  **REPORT**  **COPY**  **ENABLE**  **DELETE**

On survey end

Jump to website

Go to page (URL)

<https://www.surveylab.com/>

Total responses limit

unlimited

Quota-full (URL)

Screen-out (URL)

Type	Description
Completed	Go to page (URL) . This is a redirection link for a completed survey. When the respondent finishes the survey, he/she will be redirected to this URL.
Screen-out	Screen-out (URL) . This is redirection for the respondents that have not qualified for the survey. It is used for the screen-out questions with End survey logic. This kind of question should be placed at the beginning of your questionnaire.
Quota full	Quota full (URL) . This is a redirection link for the research project completed when a required number of responses is collected.

4. Set option Response **Key: URL** only and **SAVE** collector settings.



5. Copy survey link, activate collector, and go back to CINT

Adding the Survey URL

1. Add the parameter **/crk/[ID]/** at the end of your survey link in Access

Survey Links Setup ?

Connect your survey tool to Access/Pro. The live link is used while your project is in field.

You can have a separate link for testing or use the same live link.

Mandatory tag

[ID]	Respondent ID
------	---------------

+ Optional tags

Live link Link template Unique links

Examples
[https://www.example.com/survey.html?id=\[ID\]](https://www.example.com/survey.html?id=[ID])
[https://www.example.com/survey.html?id=\[ID\]&gender=\[S\]](https://www.example.com/survey.html?id=[ID]&gender=[S])

Link template *

Test link Link template Unique links

Link template *

Screenout and testing instructions *

Cancel
Confirm



2. To automatically add data like gender or age to the survey report use http referrer or user attributes functionality. To do it add parameter on the end of your survey link
 - a. /crk/[ID]/referrer/[S]-[YOB]-[AGE]/custom1/[S]/custom2/[YOB]/custom3/[AGE]/

b. **EXAMPLE:**

[https://www.surveylab.com/pageTag/SurveyCampaign/cld/1a74f0bf67bfe46ad22d4d9301b4b37bab47075ec/crk/\[ID\]/referrer/\[S\]-\[YOB\]-\[AGE\]/custom1/\[S\]/custom2/\[YOB\]/custom3/\[AGE\]/](https://www.surveylab.com/pageTag/SurveyCampaign/cld/1a74f0bf67bfe46ad22d4d9301b4b37bab47075ec/crk/[ID]/referrer/[S]-[YOB]-[AGE]/custom1/[S]/custom2/[YOB]/custom3/[AGE]/)

Parameter (URL)	Description
/crk/[ID]/	Response Key.
/crk/[ID]/referrer/[S]-[YOB]-[AGE]/	Response Key and referrer.
/crk/[ID]/custom1/[S]/custom2/[YOB]/custom3/[AGE]/	Response Key and user attributes.
/crk/[ID]/referrer/[S]-[YOB]-[AGE]/custom1/[S]/custom2/[YOB]/custom3/[AGE]/	Response Key, referrer and user attributes.

3. Using parameters in the survey link CINT can automatically transfer the following information into your SurveyLab survey:

Parameter	Description
[ID]	= Respondent ID
[M]	= Member ID
[P]	= Panelist ID
[PC]	= Postal code
[S]	= Gender (male = 1, female = 2)
[YOB]	= Year of birth
[AGE]	= Age